

Datasheet

KIAA0226 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009711-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant KIAA0226.

Immunogen: KIAA0226 (AAH33615.1, 1 a.a. ~ 375 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MQSILYHGLIRDQACRRQTDYWQFVKDIRWLSPHSAL
HVEKFISVHENDQSSADGASERAVAELWLQHSLQYHC
LSAQLRPLLGDRQYIRKFYTDAFLLSDAHVTAMLQCL
EAVEQNNPRLLAQIDASMFARKHESPLLVTKSQSLTAL
PSSTYTTPNSYAQHSYFGSFSSLHQSVPNNGSERRST
SFPLSGPPRKPQESRGHVSPAEDQTIQAPPVSVSALA
RDSPLTPNEMSSSTLTSPIEASWVSSQNDSPGDASEG
PEYLAIGNLDPRGRTASCQSHSSNAESSSSNLFSSSS
SQKPDSAASSLGDQEGGGESQLSSVLRSSSFSEGQT
LTVTSGAKKSHIRSHSDTSIASRGAPGNEEHRLLVSKM
TLN

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9711

Gene Symbol: KIAA0226

Gene Alias: -