

Datasheet

GPR126 polyclonal antibody

Catalog Number: PAB25978

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GPR126.

Immunogen: A synthetic peptide corresponding to 18 amino acids at internal region of human GPR126.

Host: Rabbit

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: BLAST analysis of the peptide immunogen showed no homology with other human proteins.

Form: Liquid

Purification: Immunoaffinity chromatography

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (10 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 57211

Gene Symbol: GPR126

Gene Alias: APG1, DREG, PS1TP2, STQTL1, VIGR

Gene Summary: This gene, which is upregulated in

human umbilical vein endothelial cells, encodes a G protein-coupled receptor. Variations in this gene can affect a person's stature. Multiple transcript variants encoding different proteins have been found for this gene. [provided by RefSeq]