



Certificate of Conformity

VWR LOW FORM BEAKERS, WITH RIM AND SPOUT

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (ml)	Height(mm)	O.D.(mm)	Pk
213-0292	100	60	50	1
213-0293	250	80	65	1
213-0294	500	110	80	1
213-0295	1000	130	100	1
213-0296	2000	180	120	1
213-0297	3000	180	150	1
213-0298	5000	210	180	1

Type: 18/8 (AISI 304) stainless steel, rust-proof

Complies with: ISO13485

Other Characteristics:

- A complete range of metal beakers for a wide variety of applications in the laboratory.
- For use with corrosive materials.
- Robust and durable
- Excellent resistance to heat
- Smooth surface for quick and easy cleaning

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **January 2016**